



GORE® EXCLUDER® AAA Endoprosthesis



GORE® EXCLUDER® Iliac Branch Endoprosthesis

CONFIGURATIONS

Trunk-Ipsilateral Leg Endoprosthesis

GO [®] C3 [®] DELIVERY SYSTEM CATALOGUE NUMBER	INTENDED AORTIC INNER DIAMETER (mm)	AORTIC ENDOPROSTHESIS DIAMETER (mm)	INTENDED ILIAC INNER DIAMETER (mm)	ILIAC ENDOPROSTHESIS DIAMETER (mm)	OVERALL DEVICE LENGTH (cm)	RECOMMENDED SHEATH SIZE (Fr)
RLT231212	19–21	23	10–11	12	12	16
RLT231214	19–21	23	10–11	12	14	16
RLT231216	19–21	23	10–11	12	16	16
RLT231218	19–21	23	10–11	12	18	16
RLT231412	19–21	23	12–13.5	14.5	12	16
RLT231414	19–21	23	12–13.5	14.5	14	16
RLT231416	19–21	23	12–13.5	14.5	16	16
RLT231418	19–21	23	12–13.5	14.5	18	16
RLT261212	22–23	26	10–11	12	12	16
RLT261214	22–23	26	10–11	12	14	16
RLT261216	22–23	26	10–11	12	16	16
RLT261218	22–23	26	10–11	12	18	16
RLT261412	22–23	26	12–13.5	14.5	12	16
RLT261414	22–23	26	12–13.5	14.5	14	16
RLT261416	22–23	26	12–13.5	14.5	16	16
RLT261418	22–23	26	12–13.5	14.5	18	16
RLT281212	24–26	28.5	10–11	12	12	18
RLT281214	24–26	28.5	10–11	12	14	18
RLT281216	24–26	28.5	10–11	12	16	18
RLT281218	24–26	28.5	10–11	12	18	18
RLT281412	24–26	28.5	12–13.5	14.5	12	18
RLT281414	24–26	28.5	12–13.5	14.5	14	18
RLT281416	24–26	28.5	12–13.5	14.5	16	18
RLT281418	24–26	28.5	12–13.5	14.5	18	18
RLT311413	27–29	31	12–13.5	14.5	13	18
RLT311415	27–29	31	12–13.5	14.5	15	18
RLT311417	27–29	31	12–13.5	14.5	17	18
RLT351414	30–32	35	12–13.5	14.5	14	18
RLT351416	30–32	35	12–13.5	14.5	16	18
RLT351418	30–32	35	12–13.5	14.5	18	18

Contralateral Leg Endoprosthesis

CATALOGUE NUMBER	INTENDED ILIAC INNER DIAMETER (mm)	ILIAC ENDOPROSTHESIS DIAMETER (mm)	CONTRALATERAL LEG LENGTH (cm)	RECOMMENDED SHEATH SIZE (Fr)
PLC121000	10–11	12	10	12
PLC121200	10–11	12	12	12
PLC121400	10–11	12	14	12
PLC141000	12–13.5	14.5	10	12
PLC141200	12–13.5	14.5	12	12
PLC141400	12–13.5	14.5	14	12
PLC161000*	13.5–14.5	16	9.5	12
PLC161200*	13.5–14.5	16	11.5	12
PLC161400*	13.5–14.5	16	13.5	12
PLC181000*	14.5–16.5	18	9.5	12
PLC181200*	14.5–16.5	18	11.5	12
PLC181400*	14.5–16.5	18	13.5	12
PLC201000*	16.5–18.5	20	9.5	12
PLC201200*	16.5–18.5	20	11.5	12
PLC201400*	16.5–18.5	20	13.5	12
PLC231000*	18.5–21.5	23	10	14
PLC231200*	18.5–21.5	23	12	14
PLC231400*	18.5–21.5	23	14	14
PLC271000*	21.5–25.0	27	10	15
PLC271200*	21.5–25.0	27	12	15
PLC271400*	21.5–25.0	27	14	15

Aortic Accessories

GO[®] DrySeal Flex Introducer Sheath

CATALOGUE NUMBER	SHEATH SIZE (Fr)	MINIMUM SHEATH ID (mm)	NOMINAL SHEATH OD (mm)	WORKING LENGTH (cm)
DSF1233	12	4.0	4.7	33
DSF1245	12	4.0	4.7	45
DSF1433	14	4.7	5.3	33
DSF1533	15	5.0	5.6	33
DSF1633	16	5.3	6.1	33
DSF1833	18	6.0	6.7	33

GO[®] Molding & Occlusion Balloon

CATALOGUE NUMBER	SIZE (Fr)	BALLOON DIAMETER (mm)	CATHETER LENGTH (cm)
MOB37	10	10–37	90

Iliac Extender Endoprosthesis

CATALOGUE NUMBER	INTENDED ILIAC INNER DIAMETER (mm)	ENDO PROSTHESIS DIAMETER (mm)	ENDO PROSTHESIS LENGTH (cm)	RECOMMENDED SHEATH SIZE (Fr)
PLL161007	8–9	10	7	12
PLL161207	10–11	12	7	12
PLL161407	12–13.5	14.5	7	12

Aortic Extender Endoprosthesis

CATALOGUE NUMBER	INTENDED AORTIC INNER DIAMETER (mm)	ENDO PROSTHESIS DIAMETER (mm)	ENDO PROSTHESIS LENGTH (cm)	RECOMMENDED SHEATH SIZE (Fr)
PLA230300	19–21	23	3.3	16
PLA260300	22–23	26	3.3	16
PLA280300	24–26	28.5	3.3	16
PLA320400	27–29	32	4.5	17
PLA360400	30–32	36	4.5	18

Iliac Branch Component**

CATALOGUE NUMBER	INTENDED EXTERNAL ILIAC VESSEL DIAMETER (mm)	DIAMETER ILIAC END (mm)	ILIAC BRANCH COMPONENT LENGTH (cm)	RECOMMENDED SHEATH SIZE (Fr)
CEB231010	6.5–9	10	10	16
CEB231210	10–11	12	10	16
CEB231410	12–13.5	14.5	10	16


Internal Iliac Component**

CATALOGUE NUMBER	INTENDED INTERNAL ILIAC VESSEL DIAMETER (mm)	DIAMETER ILIAC END (mm)	INTERNAL ILIAC COMPONENT LENGTH (cm)	RECOMMENDED SHEATH SIZE (Fr x cm)
HGB161007	6.5–9	10	7	12 x 45
HGB161207	10–11	12	7	12 x 45
HGB161407	12–13.5	14.5	7	12 x 45


* PLEASE NOTE: All large diameter Contralateral Leg Endoprostheses (16, 18, 20, 23, 27 mm) can be used as Iliac Extenders.

** For use in the United States, add an 'A' at the end of the catalogue number (i.e., CEB231010A).


■ **23 / 26 / 28.5 mm Aortic Extender**
(16 Fr Sheath Required)






■ **32 / 36 mm Aortic Extender**
(17 / 18 Fr Sheath Required)



■ **23 / 26 / 28.5 mm Trunk-Ipsilateral Leg Endoprosthesis**
(16 / 18 Fr Sheath Required)




■ **31 / 35 mm Trunk-Ipsilateral Leg Endoprosthesis**
(18 Fr Sheath Required)




■ **Contralateral Leg Endoprosthesis²**
(12 Fr Sheath Required)

■ **Iliac Extender¹**
(12 Fr Sheath Required)





■ **Contralateral Leg Endoprosthesis**
(12 Fr Sheath Required)



- Trunk-Ipsilateral Leg
- Contralateral Leg
- Iliac Extender
- Aortic Extender
- Iliac Branch Component
- Internal Iliac Component

¹ Traditional sizes of Iliac Extender (10, 12, 14.5 mm) are not compatible with large diameter (16, 18, 20, 23, 27 mm) Contralateral Leg Endoprostheses.
² All large diameter Contralateral Leg Endoprostheses (16, 18, 20, 23, 27 mm) can be used as Iliac Extenders.
³ The 27 mm x 10 cm Contralateral Leg Endoprosthesis does not have a radiopaque marker 6 cm above the distal end.
⁴ 12 Fr x 45 cm flexible reinforced sheath.



■ **Contralateral Leg Endoprosthesis^{2,3}**
(14 / 15 Fr Sheath Required)

■ **Internal Iliac Component⁴**
(12 Fr Sheath Required)

■ **Iliac Branch Component**
(16 Fr Sheath Required)



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)

00800.6334.4673 (Europe)

800.437.8181 (United States)

928.779.2771 (United States)

goremedical.com

Products listed may not be available in all markets.

GORE®, C3®, EXCLUDER®, and designs are trademarks of W. L. Gore & Associates.
© 2018 W. L. Gore & Associates, Inc. AX0215-EN3 SEPTEMBER 2018